


## Anti-Trpm8 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z2W7

Item Number: CM0016752



### Introduction

TRPM8 Rabbit pAb is a polyclonal antibody generated in rabbit against a recombinant protein of human TRPM8. Validated for WB and ELISA, it cross-reacts with human and mouse TRPM8 (128 kDa). Ideal for cold-sensing and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | TRPM8 Rabbit pAb                                                                                                    |
| Remarks/Alias    | TRPP8; LTRPC6; trp-p8; LTrpC-6; TRPM8                                                                               |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 79054                                                                                                               |
| GeneID           | 79054                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 880-960 of human TRPM8 (NP_076985.4). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                        |
| SWISS            | Q7Z2W7                                                                                                              |
| Protein Weight   | 128kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                             |