

Anti-Trps1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhf7

Item Number: CM0028198

Introduction

Rabbit polyclonal antibody targeting TRPS1, a zinc finger transcriptional repressor. Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for developmental biology and skeletal research applications.

[Learn More](#)

Parameter	Value
Product Name	TRPS1 Rabbit pAb
Remarks/Alias	GC79; LGCR; TRPS1
Species	Human
GeneID (Human)	7227
GeneID	7227
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-260 of human TRPS1 (NP_054831.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UHF7
Protein Weight	142kDa
Shipping	Ice bag