


Anti-Trpv1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ner1

Item Number: CM0029966



Introduction

High-quality rabbit polyclonal antibody against TRPV1 (capsaicin receptor, VR1) for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for studies on pain sensing, thermosensation, and ion channel signaling.

[Learn More](#)

Parameter	Value
Product Name	TRPV1 Rabbit pAb
Remarks/Alias	VR1; TRPV1
Species	Human
GeneID	7442 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 725-839 of human TRPV1 (NP_061197.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NER1
Protein Weight	95kDa
Shipping	Ice bag