


Anti-Trpv1 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Ner1

Item Number: CM0007184



Introduction

Rabbit monoclonal antibody targeting human TRPV1, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for pain and sensory neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	TRPV1 Rabbit mAb
Remarks/Alias	VR1; TRPV1
Species	Human
Gene ID	7442
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 725-839 of human TRPV1 (NP_542435.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8NER1
Protein Weight	95kDa
Shipping	Ice bag