


Anti-Trpv3 Polyclonal Antibody For Wb, Elisa - Q8Net8

Item Number: CM0028200



Introduction

TRPV3 Rabbit pAb targets human TRPV3 (SWISS: Q8NET8), validated for WB and ELISA with cross-reactivity to mouse. Ideal for thermosensation, ion channel, and dermatology research.

[Learn More](#)

Parameter	Value
Product Name	TRPV3 Rabbit pAb
Remarks/Alias	OLMS; VRL3; OLMS1; FNEPPK2; TRPV3
Species	Human
Gene ID	162514
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPV3 (NP_659505.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8NET8
Protein Weight	91kDa
Shipping	Ice bag