


Anti-Trpv4 Monoclonal Antibody For Wb & Elisa - Q9Hba0

Item Number: CM0007439



Introduction

TRPV4 Rabbit Monoclonal Antibody validated for WB and ELISA detects human TRPV4, a non-selective calcium channel involved in osmotic sensitivity, mechanosensation, and joint homeostasis – UniProt Q9HBA0.

[Learn More](#)

Parameter	Value
Product Name	TRPV4 Rabbit mAb
Remarks/Alias	SMAL; VRL2; BCYM3; CMT2C; SPSMA; TRP12; VROAC; HMSN2C; OTRPC4; SSQTL1; TRPV4
Species	Human
Gene ID (Human)	59341
Gene ID	59341
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 771-871 of human TRPV4 (NP_067638.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HBA0
Protein Weight	98kDa
Shipping	Ice bag