


Anti-Trpv4 Polyclonal Antibody For Wb And Elisa - Q9Hba0

Item Number: CM0027152



Introduction

Rabbit polyclonal antibody targeting human TRPV4, raised against a synthetic peptide within the C-terminal region. Suitable for Western blot and ELISA applications; validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	TRPV4 Rabbit pAb
Remarks/Alias	SMAL; VRL2; BCYM3; CMT2C; SPSMA; TRP12; VROAC; HMSN2C; OTRPC4; SSQTL1; TRPV4
Species	Human
Gene ID	59341
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 800-871 of human TRPV4 (NP_067638.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HBA0
Protein Weight	98kDa
Shipping	Ice bag