


Anti-Trrap Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y4A5

Item Number: CM0019854



Introduction

TRRAP rabbit polyclonal antibody (Q9Y4A5) validated for WB, IHC-P, ELISA. Recognizes human, rat TRRAP (~438 kDa). Essential for chromatin remodeling and MYC transcription research.

[Learn More](#)

Parameter	Value
Product Name	TRRAP Rabbit pAb
Remarks/Alias	Tra1; TR-AP; DEDDFA; DFNA75; PAF400; STAF40; PAF350/400; TRRAP
Species	Human
GeneID (Human)	8295
GeneID	8295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3660-3859 of TRRAP (NP_001231509.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9Y4A5
Protein Weight	438kDa
Shipping	Ice bag