


Anti-Trypsin Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P07477

Item Number: CM0011836



Introduction

Rabbit monoclonal antibody against human trypsin (PRSS1), validated for immunofluorescence, immunohistochemistry, and ELISA. Recognizes human trypsin; ideal for digestive enzyme and protease research. Immunogen: recombinant fusion protein (aa 24-247).

[Learn More](#)

Parameter	Value
Product Name	Trypsin Rabbit mAb
Remarks/Alias	TRP1; TRY1; TRY4; TRYP1
Species	Human
GeneID (Human)	5644
GeneID	5644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-247 of human PRSS1 (NP_002760.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P07477
Protein Weight	26 kDa
Shipping	Ice bag