

Anti-Tryptophanyl-Trna Synthetase 1 Monoclonal Antibody For Wb & Elisa - P23381

Item Number: CM0011247



Introduction

Rabbit monoclonal antibody targeting human Tryptophanyl-tRNA synthetase 1 (WARS1). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for translation and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Tryptophanyl-tRNA synthetase 1 Rabbit mAb
Remarks/Alias	HMN9; WARS; IFI53; IFP53; GAMMA-2; Tryptophanyl-tRNA synthetase 1
Species	Human
Gene ID (Human)	7453
Gene ID	7453
Immunogen	Recombinant fusion protein containing amino acids 3-98 of human Tryptophanyl-tRNA synthetase 1 (P23381)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23381
Protein Weight	53kDa
Shipping	Ice bag