

Anti-Tryptophanyl-Trna Synthetase 1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P23381

Item Number: CM0007232



Introduction

High-quality rabbit monoclonal antibody targeting human Tryptophanyl-tRNA synthetase 1 (WARS1). Validated for WB, IHC-P, ELISA. Ideal for protein synthesis, angiogenesis, and interferon signaling research.

[Learn More](#)

Parameter	Value
Product Name	Tryptophanyl-tRNA synthetase 1 Rabbit mAb
Remarks/Alias	HMN9; WARS; IFI53; IFP53; GAMMA-2; Tryptophanyl-tRNA synthetase 1
Species	Human
GeneID (Human)	7453
GeneID	7453
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 372-471 of human Tryptophanyl-tRNA synthetase 1 (P23381)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P23381
Protein Weight	53kDa
Shipping	Ice bag