


Anti-Tsc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q92574

Item Number: CM0014066



Introduction

Rabbit polyclonal anti-TSC1 (Hamartin) antibody for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for mTOR pathway research and tumor suppressor analysis.

[Learn More](#)

Parameter	Value
Product Name	TSC1 Rabbit pAb
Remarks/Alias	LAM; TSC; TSC1
Species	Human
GeneID (human)	7248
GeneID	7248
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 600-700 of human TSC1 (NP_000359.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92574
Protein Weight	130kDa
Shipping	Ice bag