


Anti-Tsc2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49815

Item Number: CM0015183



Introduction

High-quality rabbit polyclonal antibody targeting human TSC2 (tuberin), validated for WB, IHC-P, IF/ICC, and ELISA. Detects TSC2 in human, mouse, and rat samples, suitable for mTORC1 signaling and tuberous sclerosis research.

[Learn More](#)

Parameter	Value
Product Name	TSC2 Rabbit pAb
Remarks/Alias	LAM; TSC4; PPP1R160; TSC2
Species	Human
GeneID (Human)	7249
GeneID	7249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1160-1255 of human TSC2 (NP_000539.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49815
Protein Weight	201kDa
Shipping	Ice bag