


Anti-Tsc2 (Tuberin) Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P49815

Item Number: CM0003082



Introduction

Rabbit monoclonal antibody targeting human TSC2 (tuberin). Validated for Western blot, immunoprecipitation, and ELISA. Reacts with human. MW ~201 kDa. Suitable for mTOR/TSC signaling research.

[Learn More](#)

Parameter	Value
Product Name	TSC2 Rabbit mAb
Remarks/Alias	LAM; TSC4; PPP1R160; TSC2
Species	Human
GeneID (human)	7249
GeneID	7249
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1700-1807 of human Tuberin (P49815).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P49815
Protein Weight	201kDa
Shipping	Ice bag