

Anti-Tsc22D3 Polyclonal Antibody For If/Icc, Elisa - Q99576

Item Number: CM0024611

Introduction

Rabbit polyclonal antibody targeting human TSC22D3 (GILZ), suitable for IF/ICC and ELISA. Cross-reacts with rat. Ideal for apoptosis, inflammation, and glucocorticoid signaling research.

[Learn More](#)

Parameter	Value
Product Name	TSC22D3 Rabbit pAb
Remarks/Alias	DIP; GILZ; DSIPI; TSC-22R; TSC22D3
Species	Human
GeneID (human)	1831
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TSC22D3 (NP_004080.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Rat
SWISS	Q99576
Protein Weight	15kDa
Shipping	Ice bag