

## Anti-Tsc22D3 Polyclonal Antibody For Wb, Elisa - Q99576

Item Number: CM0029722

### Introduction

High-quality rabbit polyclonal antibody targeting human TSC22D3 (GILZ, DSIPI). Tested in WB and ELISA. Ideal for studying apoptosis, inflammation, and glucocorticoid signaling. Reliable for IVD raw material development.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | TSC22D3 Rabbit pAb                                                                                                  |
| Remarks/Alias    | DIP; GILZ; DSIPI; TSC-22R; TSC22D3                                                                                  |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 1831                                                                                                                |
| GeneID           | 1831                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 58-134 of human TSC22D3 (NP_004080.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | Q99576                                                                                                              |
| Protein Weight   | 15kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |