

Anti-Tsen2 Polyclonal Antibody For Wb, Elisa - Q8Nce0

Item Number: CM0029700



Introduction

High-quality rabbit polyclonal antibody against human TSEN2 (tRNA-splicing endonuclease subunit Sen2). Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	TSEN2 Rabbit pAb
Remarks/Alias	SEN2; PCH2B; SEN2L; TSEN2
Species	Human
Gene ID (Human)	80746
Gene ID	80746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TSEN2 (NP_079541.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NCE0
Protein Weight	53kDa
Shipping	Ice bag