

Anti-Tsen54 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z6J9

Item Number: CM0017065



Introduction

Rabbit polyclonal antibody targeting TSEN54, a tRNA-splicing endonuclease subunit involved in RNA processing and neurodegeneration. Validated for WB and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	TSEN54 Rabbit pAb
Remarks/Alias	PCH4; PCH5; PCH2A; sen54; SEN54L; TSEN54
Species	Human
GeneID (Human)	283989
GeneID	283989
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TSEN54 (NP_997229.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z6J9
Protein Weight	59kDa
Shipping	Ice bag