


Anti-Tsfm Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P43897

Item Number: CM0012838



Introduction

Rabbit polyclonal antibody against human TSFM (mitochondrial elongation factor Ts). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial translation and disease research.

[Learn More](#)

Parameter	Value
Product Name	TSFM Rabbit pAb
Remarks/Alias	EFTS; EFTSMT; TSFM
Species	Human
Gene ID (Human)	10102
Gene ID	10102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-346 of human TSFM (NP_001166167.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43897
Protein Weight	35kDa
Shipping	Ice bag