

Anti-Tsg101/Vps23 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99816

Item Number: CM0024057



Introduction

Rabbit polyclonal antibody against human TSG101/VPS23, validated for WB, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Ideal for ESCRT-I complex, exosome, and viral budding research.

[Learn More](#)

Parameter	Value
Product Name	TSG101/VPS23 Rabbit pAb
Remarks/Alias	TSG10, VPS23, TSG101/VPS23
Species	Human
GeneID (Human)	7251
GeneID	7251
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 291-390 of human TSG101/VPS23 (NP_006283.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99816
Protein Weight	44kDa
Shipping	Ice bag