

Anti-Tsg101/Vps23 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99816

Item Number: CM0019779



Introduction

High-quality rabbit polyclonal antibody against TSG101/VPS23 (ESCRT-I subunit). Validated for WB, IF/ICC, and ELISA; reacts with human, mouse, and rat. Ideal for endosomal trafficking, viral budding, and exosome research.

[Learn More](#)

Parameter	Value
Product Name	TSG101/VPS23 Rabbit pAb
Remarks/Alias	TSG10; VPS23; TSG101/VPS23
Species	Human
GeneID (Human)	7251
GeneID	7251
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 291-390 of human TSG101/VPS23 (NP_006283.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99816
Protein Weight	44kDa
Shipping	Ice bag