

Anti-Tshb Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P01222

Item Number: CM0026041



Introduction

TSHB Rabbit Polyclonal Antibody validated for IHC-P, IF/ICC, ELISA. Detects human and mouse TSH-beta, ~16 kDa. Ideal for thyroid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TSHB Rabbit pAb
Remarks/Alias	TSH-B; TSH-BETA; TSHB
Species	Human
GeneID (Human)	7252
GeneID	7252
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-138 of human TSHB (NP_000540.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01222
Protein Weight	16kDa
Shipping	Ice bag