

Anti-Tshb Rabbit Polyclonal Antibody For Wb, Elisa - P01222

Item Number: CM0015625



Introduction

Rabbit polyclonal antibody targeting human TSHB (TSH-beta), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for thyroid hormone and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TSHB Rabbit pAb
Remarks/Alias	TSH-B; TSH-BETA; TSHB
Species	Human
Gene ID (Human)	7252
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human TSHB (NP_000540.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01222
Protein Weight	16kDa
Shipping	Ice bag