

Anti-Tshr Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P16473

Item Number: CM0028019



Introduction

CamelBio's TSHR Rabbit Polyclonal Antibody detects human, mouse, and rat TSHR for western blot, IHC-P, and ELISA. Ideal for thyroid hormone signaling and receptor research.

[Learn More](#)

Parameter	Value
Product Name	TSHR Rabbit pAb
Remarks/Alias	LGR3; CHNG1; hTSHR-I; TSHR
Species	Human
Gene ID (Human)	7253
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-253 of human TSHR (NP_001018046.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16473
Protein Weight	87kDa
Shipping	Ice bag