


Anti-Tshz3 Rabbit Polyclonal Antibody For Wb, Elisa - Q63Hk5

Item Number: CM0023087



Introduction

High-quality TSHZ3 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human and mouse TSHZ3 protein, a transcriptional regulator crucial for neurodevelopment and ureter smooth muscle formation. Suitable for neuroscience and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	TSHZ3 Rabbit pAb
Remarks/Alias	TSH3; ZNF537; TSHZ3
Species	Human
GeneID (Human)	57616
GeneID	57616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 791-890 of human TSHZ3 (NP_065907.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q63HK5
Protein Weight	119kDa
Shipping	Ice bag