


Anti-Tslp Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - Q969D9

Item Number: CM0007391



Introduction

Recombinant Rabbit monoclonal TSLP antibody (SWISS: Q969D9) for WB, IHC, IF, ELISA. Cross-reacts with human and mouse. Ideal for dendritic cell maturation and allergic inflammation research.

[Learn More](#)

Parameter	Value
Product Name	TSLP Rabbit mAb
Remarks/Alias	TSLP
Species	Human
GeneID (Human)	85480
GeneID	85480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human TSLP(NP_149024.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q969D9
Protein Weight	18kDa
Shipping	Ice bag