


Anti-Tslp Rabbit Polyclonal Antibody For Wb, If/Icc, And Elisa - Q969D9

Item Number: CM0018588



Introduction

Rabbit polyclonal antibody against human TSLP, suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for investigating TSLP-mediated dendritic cell maturation and allergic inflammation.

[Learn More](#)

Parameter	Value
Product Name	TSLP Rabbit pAb
Remarks/Alias	TSLP
Species	Human
GeneID	85480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human TSLP (NP_149024.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969D9
Protein Weight	18kDa
Shipping	Ice bag