


Anti-Tsnax Rabbit Polyclonal Antibody For Wb, Elisa - Q99598

Item Number: CM0016324



Introduction

Rabbit polyclonal antibody against human TSNAX (TRAX, C3PO), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Supports RNA silencing complex and spermatogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	TSNAX Rabbit pAb
Remarks/Alias	C3PO; TRAX; TSNAX
Species	Human
GeneID	7257
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human TSNAX (NP_005990.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99598
Protein Weight	33kDa
Shipping	Ice bag