

Anti-Tspan31 Polyclonal Antibody For Wb, Elisa - Q12999

Item Number: CM0015655

CamelBio is a leading provider of high-quality antibodies and reagents for research and clinical applications. Our products are sourced from reputable suppliers and are rigorously tested to ensure purity and performance. We offer a wide range of products, including monoclonal and polyclonal antibodies, enzymes, and various reagents, all designed to meet the needs of our customers. Our commitment to quality and customer satisfaction is reflected in our extensive product catalog and our responsive customer support team.

Introduction

Rabbit polyclonal antibody targeting human TSPAN31 (SAS), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for membrane protein research, cancer studies, and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	TSPAN31 Rabbit pAb
Remarks/Alias	SAS; TSPAN31
Species	Human
GeneID (Human)	6302
GeneID	6302
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 95-175 of human TSPAN31 (NP_005972.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12999
Protein Weight	23kDa
Shipping	Ice bag