

Anti-Tspan33 Polyclonal Antibody For Wb, Elisa - Q86Uf1

Item Number: CM0029712

Introduction

TSPAN33 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human, mouse, and rat TSPAN33. Synthetic peptide immunogen targeting the C-terminus. Ideal for erythropoiesis and ADAM10 interaction studies.

[Learn More](#)

Parameter	Value
Product Name	TSPAN33 Rabbit pAb
Remarks/Alias	PEN; PEN.; TSPAN-33; TSPAN33
Species	Human
Gene ID (Human)	340348
Gene ID	340348
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human TSPAN33 (NP_848657.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UF1
Protein Weight	32kDa
Shipping	Ice bag