

Anti-Tspan4/Nag-2 Rabbit Polyclonal Antibody For Wb, Elisa - O14817

Item Number: CM0017471

Introduction

High-quality rabbit polyclonal antibody targeting human TSPAN4/NAG-2 (also known as NAG2, TM4SF7). Validated for WB and ELISA, with cross-reactivity in mouse and rat. Ideal for studying tetraspanin-mediated membrane microdomains and migrasome biology.

[Learn More](#)

Parameter	Specification
Product Name	TSPAN4/NAG-2 Rabbit pAb
Remarks/Alias	NAG2; NAG-2; TM4SF7; TSPAN-4; TETRASPAN; TSPAN4/NAG-2
Species	Human
GeneID (Human)	7106
GeneID	7106
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human TSPAN4/NAG-2 (NP_003262.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14817
Protein Weight	26kDa
Shipping	Ice bag