

Anti-Tspan7 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P41732

Item Number: CM0026188

Introduction

TSPAN7 Rabbit pAb (Polyclonal Antibodies) validated for WB, IF-P, ELISA. Detects human TSPAN7 (28kDa) with cross-reactivity to mouse and rat. Ideal for cell proliferation and disease research.

[Learn More](#)

Parameter	Value
Product Name	TSPAN7 Rabbit pAb
Remarks/Alias	A15; MXS1; CD231; MRX58; CCG-B7; TM4SF2; XLID58; TALLA-1; TM4SF2b; DXS1692E; TSPAN7
Species	Human
GeneID	7102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-240 of human TSPAN7 (NP_004606.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P41732
Protein Weight	28kDa
Shipping	Ice bag