

Anti-Tspo Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P30536

Item Number: CM0022287

Introduction

Rabbit polyclonal antibody against human TSPO (PBR), validated for Western blot, IF/ICC, and ELISA, cross-reacts with mouse. Detects the 19 kDa mitochondrial translocator protein involved in cholesterol transport and steroid synthesis.

[Learn More](#)

Parameter	Value
Product Name	PBR/TSPO Rabbit pAb
Remarks/Alias	DBI; IBP; MBR; PBR; PBS; BPBS; BZRP; PKBS; PTBR; mDRC; pk18; TSPO1; PBR/TSPO
Species	Human
GeneID (Human)	706
GeneID	706
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-50 of human PBR/PBR/TSPO (NP_000705.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P30536
Protein Weight	19kDa
Shipping	Ice bag