

Anti-Tspyl2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H2G4

Item Number: CM0013036

Introduction

Rabbit polyclonal antibody targeting human TSPYL2 (Testis-specific Y-encoded-like protein 2), validated for Western blotting and ELISA. Cross-reacts with mouse and rat TSPYL2. Suitable for investigating transcriptional regulation, chromatin dynamics, and cell cycle control.

[Learn More](#)

Parameter	Value
Product Name	TSPYL2 Rabbit pAb
Remarks/Alias	CDA1; CTCL; NP79; TSPX; CINAP; DENTT; SE204; HRIHFB2216; TSPYL2
Species	Human
Gene ID (Human)	64061
Gene ID	64061
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 240-460 of human TSPYL2 (NP_071400.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2G4
Protein Weight	79kDa
Shipping	Ice bag