

Anti-Tssk3 Polyclonal Antibody For Wb, Elisa - Q96Pn8

Item Number: CM0031086

Introduction

TSSK3 Rabbit pAb validated for WB and ELISA with cross-reactivity to mouse and rat. Targets testis-specific serine/threonine kinase 3, essential for spermatid development and male fertility. Ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	TSSK3 Rabbit pAb
Remarks/Alias	TSK3; SPOGA3; STK22C; STK22D; TSSK3
Species	Human
GeneID (Human)	81629
GeneID	81629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human TSSK3 (NP_443073.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96PN8
Protein Weight	30kDa
Shipping	Ice bag