

Anti-Tssk4 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Sa08

Item Number: CM0030735

Introduction

Rabbit polyclonal antibody targeting human TSSK4, validated for WB and ELISA; cross-reacts with mouse. The ~37 kDa testis-specific serine/threonine kinase is implicated in spermatogenesis and Cre/Creb signaling. Research use.

[Learn More](#)

Parameter	Value
Product Name	TSSK4 Rabbit pAb
Remarks/Alias	TSK4; TSK-4; TSSK5; STK22E; C14orf20; TSSK4
Species	Human
Gene ID (Human)	283629
Gene ID	283629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 79-328 of human TSSK4 (NP_777604.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q6SA08
Protein Weight	37kDa
Shipping	Ice bag