


Anti-Ttbk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Iq55

Item Number: CM0026825



Introduction

Rabbit polyclonal antibody against human TTBK2, validated for WB, IHC-P, and ELISA. Cross-reacts with rat. Target is tau-tubulin kinase 2, a regulator of ciliogenesis and tau phosphorylation. Suitable for neurodegeneration and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	TTBK2 Rabbit pAb
Remarks/Alias	TTBK; SCA11; TTBK2
Species	Human
Gene ID (Human)	146057
Gene ID	146057
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1185-1244 of human TTBK2 (NP_775771.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6IQ55
Protein Weight	137kDa
Shipping	Ice bag