


Anti-Ttc11/Fis1 Polyclonal Antibody For Wb, Ip, If/Icc, Elisa - Q9Y3D6

Item Number: CM0026400



Introduction

Rabbit polyclonal antibody targeting TTC11/FIS1, suitable for WB, IP, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial dynamics, peroxisomal fission, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	TTC11/FIS1 Rabbit pAb
Remarks/Alias	TTC11; CGI-135; TTC11/FIS1
Species	Human
GeneID (Human)	51024
GeneID	51024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human TTC11/FIS1 (NP_057152.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3D6
Protein Weight	17 kDa
Shipping	Ice bag