

Anti-Ttc3 Rabbit Polyclonal Antibody For Wb, Elisa - P53804

Item Number: CM0012958

Introduction

Rabbit polyclonal antibody targeting TTC3 E3 ubiquitin-protein ligase. Validated for WB and ELISA; cross-reacts with human, mouse, rat. Suitable for ubiquitination, Akt signaling, and neuronal differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	TTC3 Rabbit pAb
Remarks/Alias	DCRR1; RNF105; TPRDIII; TTC3
Species	Human
GeneID (Human)	7267
GeneID	7267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1700 of human TTC3 (NP_003307.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53804
Protein Weight	230kDa
Shipping	Ice bag