


Anti-Ttc5 Polyclonal Antibody For Wb, Elisa - Q8N0Z6

Item Number: CM0014097



Introduction

Rabbit polyclonal antibody targeting human TTC5 (Strap), validated for Western blot and ELISA. Cross-reacts with human, mouse, rat. Ideal for chromatin regulation and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	TTC5 Rabbit pAb
Remarks/Alias	Strap; NEDCAFD; TTC5
Species	Human
GeneID (Human)	91875
GeneID	91875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TTC5 (NP_612385.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N0Z6
Protein Weight	49kDa
Shipping	Ice bag