

Anti-Ttc7A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ult0

Item Number: CM0023069

Introduction

Rabbit polyclonal antibody against human TTC7A, suitable for Western blot and ELISA applications. Cross-reacts with rat. Ideal for research on PI4K plasma membrane localization and phosphatidylinositol 4-phosphate synthesis.

[Learn More](#)

Parameter	Description
Product Name	TTC7A Rabbit pAb
Remarks/Alias	TTC7; GIDID; MINAT; TTC7A
Species	Human
GeneID	57217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TTC7A (NP_001275880.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9ULT0
Protein Weight	96kDa
Shipping	Ice bag