

Anti-Ttf-1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P43699

Item Number: CM0009927

Introduction

Anti-TTF-1 Rabbit mAb validated for WB, IF-P, IHC-P, ELISA with human, mouse, rat reactivity. Targets NKX2-1, a thyroid transcription factor crucial for lung and thyroid development. Ideal for cancer and developmental research.

[Learn More](#)

Parameter	Value
Product Name	TTF-1 Rabbit mAb
Remarks/Alias	BCH; BHC; NK-2; TEBP; TTF1; NKX2A; NMTCl; T/EBP; TITF1; TTF-1; NKX2.1; NKX2-1
Species	Human
GeneID (Human)	7080
GeneID	7080
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NKX2-1 (P43699).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43699
Protein Weight	39kDa
Shipping	Ice bag