

Anti-Ttf1 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q15361

Item Number: CM0024487

Introduction

High-quality rabbit polyclonal antibody targeting human TTF1, a nucleolar transcription termination factor. Validated for IF/ICC and ELISA. Ideal for studies on rDNA transcription regulation and nucleolar function.

[Learn More](#)

Parameter	Value
Product Name	TTF1 Rabbit pAb
Remarks/Alias	TTF-1; TTF-I; TTF1
Species	Human
GeneID	7270
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 806-905 of human TTF1 (NP_031370.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15361
Protein Weight	103kDa
Shipping	Ice bag