

Anti-Ttf2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uny4

Item Number: CM0029835

Introduction

Rabbit polyclonal antibody against human TTF2, a transcription termination factor involved in RNA polymerase II removal. Validated for WB, IF/ICC, and ELISA, supporting research in gene regulation and splicing.

[Learn More](#)

Parameter	Value
Product Name	TTF2 Rabbit pAb
Remarks/Alias	F2; HuF2; ZGRF6; TTF2
Species	Human
GeneID (Human)	8458
GeneID	8458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TTF2 (NP_003585.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UNY4
Protein Weight	130kDa
Shipping	Ice bag