

Anti-Ttk Polyclonal Antibody For Wb, If/Icc, Elisa - P33981

Item Number: CM0024033



Introduction

Rabbit polyclonal TTK (MPS1) antibody for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat TTK; targets mitotic spindle assembly checkpoint kinase. Ideal for cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	TTK Rabbit pAb
Remarks/Alias	ESK; PYT; CT96; MPH1; MPS1; MPS1L1; TTK
Species	Human
Gene ID (Human)	7272
Gene ID	7272
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 536-789 of human TTK (NP_003309.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33981
Protein Weight	97kDa
Shipping	Ice bag