

Anti-Ttll12 Polyclonal Antibody For Wb And Elisa - Q14166

Item Number: CM0023234

Introduction

Rabbit polyclonal antibody targeting human TTLL12, validated for Western blot and ELISA. Cross-reacts with mouse and rat. TTLL12 is a regulator of tubulin post-translational modifications, mitosis, chromosome stability, and innate antiviral immunity.

[Learn More](#)

Parameter	Value
Product Name	TTLL12 Rabbit pAb
Remarks/Alias	dj526l14.2; TTLL12
Species	Human
Gene ID	23170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 485-644 of human TTLL12 (NP_055955.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14166
Protein Weight	74kDa
Shipping	Ice bag