


Anti-Ttr Rabbit Monoclonal Antibody For Wb, Elisa - P02766

Item Number: CM0011584



Introduction

Anti-TTR Rabbit Monoclonal Antibody validated for WB and ELISA. Targets human transthyretin (P02766), a thyroid hormone transporter. Ideal for studying amyloidosis and hormone transport mechanisms.

[Learn More](#)

Parameter	Value
Product Name	TTR Rabbit mAb
Remarks/Alias	CTS; TTN; ATTR; CTS1; PALB; TBPA; HEL111; HsT2651; TTR
Species	Human
GeneID (Human)	7276
GeneID	7276
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 68-147 of human TTR (P02766).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02766
Protein Weight	16kDa
Shipping	Ice bag