

Anti-Tubb8 Polyclonal Antibody For Wb, If/Icc, Elisa - Q3Zcm7

Item Number: CM0014992

Introduction

High-quality rabbit polyclonal antibody against human TUBB8, validated for WB, IF/ICC, and ELISA. Ideal for studying meiotic spindle assembly and oocyte maturation in reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	TUBB8 Rabbit pAb
Remarks/Alias	OOMD; OOMD2; OZEMA2; bA631M21.2; TUBB8
Species	Human
Gene ID (Human)	347688
Gene ID	347688
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 344-444 of human TUBB8 (NP_817124.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q3ZCM7
Protein Weight	50kDa
Shipping	Ice bag