

Anti-Tubgcp3 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Cw5

Item Number: CM0029740



Introduction

Rabbit polyclonal antibody against gamma-tubulin complex component 3 (TUBGCP3/GCP3), validated for Western blot (WB) and ELISA. Reacts with human, mouse, rat. Suitable for centrosome and microtubule nucleation studies.

[Learn More](#)

Parameter	Value
Product Name	TUBGCP3 Rabbit pAb
Remarks/Alias	104p; ALP6; GCP3; Spc98; SPBC98; Spc98p; Grip104; TUBGCP3
Species	Human
Gene ID (Human)	10426
Gene ID	10426
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TUBGCP3 (NP_001273207.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96CW5
Protein Weight	104kDa
Shipping	Ice bag