

Anti-Tubgcp5 Polyclonal Antibody For Wb, Elisa - Q96Rt8

Item Number: CM0020646



Introduction

TUBGCP5 rabbit polyclonal antibody (118 kDa) for WB and ELISA applications. Detects human and rat TUBGCP5, a key component of the gamma-tubulin ring complex involved in microtubule nucleation. Suitable for centrosome biology research.

[Learn More](#)

Parameter	Value
Product Name	TUBGCP5 Rabbit pAb
Remarks/Alias	GCP5; TUBGCP5
Species	Human
GeneID (human)	114791
GeneID	114791
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human TUBGCP5 (NP_443135.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q96RT8
Protein Weight	118kDa
Shipping	Ice bag